

2'-O-Methylcytidine (Cm) Rabbit Monoclonal Antibody For Dot Blot And Elisa

Item Number: CM0000038

Introduction

Rabbit monoclonal antibody specific to 2'-O-methylcytidine (Cm), a modified RNA nucleoside. Suitable for Dot Blot (DB) and ELISA applications with species-independent cross reactivity. Ideal for epitranscriptomics research.

[Learn More](#)

Parameter	Value
Product Name	2'-O-methylcytidine/ Cm Rabbit mAb
Remarks/Alias	Cm; 2'-O-methylcytidine; 2'-O-methylcytidine/ Cm
Species	Human
Immunogen	Chemical compounds conjugated to Cm
Source	Rabbit
Category	Other Modified Antibodies
Application	Dot Blot (DB), ELISA
Cross Reactivity	Species independent
Shipping	Ice bag