


Aldoc Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09972

Item Number: CM0013906



Introduction

Rabbit polyclonal antibody against human ALDOC (fructose-bisphosphate aldolase C). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. For glycolysis and gluconeogenesis research.

[Learn More](#)

Parameter	Value
Product Name	ALDOC Rabbit pAb
Remarks/Alias	ALDC; ALDOC
Species	Human
GeneID (human)	230
GeneID	230
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human ALDOC (NP_005156.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09972
Protein Weight	39kDa
Shipping	Ice bag