

Alpk1 Rabbit Mab For Wb, Ihc-P, Elisa - E9Qnx4

Item Number: CM0007373

Introduction

High-quality ALPK1 rabbit monoclonal antibody validated for Western blot, IHC-P, and ELISA. Reacts with mouse and rat ALPK1 (139 kDa), ideal for kinase research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	ALPK1 Rabbit mAb
Remarks/Alias	Alpha-kinase 1; ALPK1
Species	Mouse
GeneID (Human)	80216
GeneID	71481
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-478 of mouse ALPK1 (NP_082084.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	E9QNX4
Protein Weight	139kDa
Shipping	Ice bag