

Anti-4F2Hc/Cd98/Slc3A2 Monoclonal Antibody For Wb, Elisa - P10852

Item Number: CM0007371



Introduction

Recombinant monoclonal antibody against mouse SLC3A2 (SWISS P10852), validated for WB and ELISA. Cross-reacts with mouse and rat. Useful for studying amino acid transporter function and integrin-related signaling.

[Learn More](#)

Parameter	Value
Product Name	4F2hc/CD98/SLC3A2 Rabbit mAb
Remarks/Alias	4F2; Cd98; Ly10; Mdu1; 4F2HC; Ly-10; NACAE; Ly-m10; Mgp-2hc; 4F2hc/CD98/SLC3A2
Species	Mouse
GeneID (Human)	6520
GeneID	17254
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 105-526 of mouse 4F2hc/CD98/SLC3A2 (P10852).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P10852
Protein Weight	58KDa
Shipping	Ice bag