


Anti-A-Actinin-4 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - O43707

Item Number: CM0008977



Introduction

Anti- α -Actinin-4 rabbit monoclonal antibody validated for WB, IF/ICC, IP, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying actin cytoskeleton, cell junctions, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	α -Actinin-4 Rabbit mAb
Remarks/Alias	FSGS; FSGS1; ACTININ-4; α -Actinin-4
Species	Human
GeneID (Human)	81
GeneID	81
Immunogen	Synthetic Peptide corresponding to amino acids 1-100 of human α -Actinin-4 (O43707)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43707
Protein Weight	105kDa
Shipping	Ice bag