

Anti-Abca8 Polyclonal Antibody For Wb, Elisa - 094911

Item Number: CM0021169

Introduction

Rabbit polyclonal antibody targeting ABCA8, validated for Western blot and ELISA in human and rat. This ATP-binding cassette transporter imports organic anions and mediates cholesterol efflux, making the antibody suitable for drug transport and lipid metabolism research.

[Learn More](#)

Parameter	Specification
Product Name	ABCA8 Rabbit pAb
Remarks/Alias	ABCA8
Species	Human
GeneID (human)	10351
GeneID	10351
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 425-495 of human ABCA8 (NP_009099.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O94911
Protein Weight	184kDa
Shipping	Ice bag