

Anti-Abcg1 Rabbit Monoclonal Antibody For Wb, Elisa - P45844

Item Number: CM0003002



Introduction

ABCG1 rabbit monoclonal antibody validated for WB and ELISA, recognizing human, mouse, and rat ABCG1. Designed for studies of cholesterol efflux and lipid transport in macrophages and other tissues.

[Learn More](#)

Parameter	Value
Product Name	ABCG1 Rabbit mAb
Remarks/Alias	ABC8; WHITE1; ABCG1
Species	Human
GeneID (Human)	9619
GeneID	9619
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45844
Protein Weight	76kDa
Shipping	Ice bag