

Anti-Abi4 Rabbit Polyclonal Antibody For Wb And Elisa - A0Mes8

Item Number: CM0019171

Introduction

Rabbit polyclonal antibody targeting Arabidopsis thaliana ABI4, an ethylene-responsive transcription factor involved in ABA signaling, stress responses, and development. Validated for WB and ELISA, ideal for plant biology research.

[Learn More](#)

Parameter	Value
Product Name	ABI4 Rabbit pAb
Remarks/Alias	ABA INSENSITIVE 4; ABSCISIC ACID-INSENSITIVE PROTEIN 4; ATABI4; GIN6; GLUCOSE INSENSITIVE 6; IMPAIRED SUCROSE INDUCTION 3; ISI3; SALOBRENO 5; SAN5; SIS5; SUCROSE UNCOUPLED 6; SUGAR-INSENSITIVE 5; SUN6; T7M7.16; ABI4
Species	Arabidopsis thaliana
GeneID (human)	64116
GeneID	818614
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of Arabidopsis thaliana ABI4 (NP_181551.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	A0MES8
Protein Weight	36kDa
Shipping	Ice bag