

Anti-Acadsb Rabbit Polyclonal Antibody For Ihc-P, Elisa - P45954

Item Number: CM0016733



Introduction

ACADSB Rabbit Polyclonal Antibody validated for IHC-P and ELISA; detects human, mouse, rat ACADSB. Ideal for studying branched-chain fatty acid oxidation and L-isoleucine metabolism. SWISS: P45954.

[Learn More](#)

Parameter	Value
Product Name	ACADSB Rabbit pAb
Remarks/Alias	ACAD7; SBCAD; 2-MEBCAD; ACADSB
Species	Human
Gene ID (Human)	36
Gene ID	36
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 133-432 of human ACADSB (NP_001600.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45954
Protein Weight	47kDa
Shipping	Ice bag