

Anti-Acc1 (Ko Validated) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q5Swu9

Item Number: CM0031090

Introduction

Rabbit polyclonal antibody targeting ACC1, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Essential for fatty acid synthesis research; supplied by CamelBio as a one-stop IVD raw material partner.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ACC1 Rabbit pAb
Remarks/Alias	Acac; Acc1; Gm738; A530025K05Rik
Species	Mouse
GeneID (human)	31
GeneID	107476
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 2246-2345 of human ACC1 (NP_579938.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5SWU9
Protein Weight	265 kDa/270 kDa
Shipping	Ice bag