

Anti-Acsl3 Polyclonal Antibody For Wb, Ihc-P, Elisa - 095573

Item Number: CM0016890

Introduction

ACSL3 Rabbit pAb validated for WB, IHC-P, ELISA; detects endogenous ACSL3 in human, mouse, rat. Long-chain-fatty-acid-CoA ligase 3, key in lipid metabolism and VLDL formation.

[Learn More](#)

Parameter	Value
Product Name	ACSL3 Rabbit pAb
Remarks/Alias	ACS3; FAFL3; LACS3; LACS 3; PRO2194; ACSL3
Species	Human
Gene ID (Human)	2181
Gene ID	2181
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 621-720 of human ACSL3 (NP_004448.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095573
Protein Weight	80kDa
Shipping	Ice bag