

## Anti-AcsL5 Rabbit Polyclonal Antibody For Wb, Elisa - Q9ULc5

Item Number: CM0013341

### Introduction

High-quality ACSL5 Rabbit Polyclonal Antibody validated for WB and ELISA with cross-reactivity in human, mouse, and rat. Ideal for studying fatty acid metabolism and lipid synthesis.

[Learn More](#)

| Parameter        | Value                                                    |
|------------------|----------------------------------------------------------|
| Product Name     | ACSL5 Rabbit pAb                                         |
| Remarks/Alias    | ACS2; ACS5; FACL5; ACSL5                                 |
| Species          | Human                                                    |
| GeneID (Human)   | 51703                                                    |
| GeneID           | 51703                                                    |
| Immunogen        | Recombinant protein (amino acids 500-739 of human ACSL5) |
| Source           | Rabbit                                                   |
| Category         | Polyclonal Antibodies                                    |
| Application      | WB, ELISA                                                |
| Cross Reactivity | Human, Mouse, Rat                                        |
| SWISS            | Q9ULC5                                                   |
| Protein Weight   | 76kDa                                                    |
| Shipping         | Ice bag                                                  |