

Anti-Adam17 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P78536

Item Number: CM0014788



Introduction

Rabbit polyclonal antibody against human ADAM17, validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for studying TNF-alpha shedding, Notch signaling, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	ADAM17 Rabbit pAb
Remarks/Alias	CSVP; TACE; NISBD; ADAM18; CD156B; NISBD1; ADAM17
Species	Human
GeneID (human)	6868
GeneID	6868
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 700-824 of human ADAM17 (NP_003174.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78536
Protein Weight	93kDa
Shipping	Ice bag