


Anti-Adgrb3 Rabbit Polyclonal Antibody For Wb, Elisa - O60242

Item Number: CM0017197



Introduction

Anti-ADGRB3 rabbit polyclonal antibody validated for WB and ELISA. Targets human Adhesion G protein-coupled receptor B3 (BAI3), a regulator of synaptogenesis and myoblast fusion. For IVD research.

[Learn More](#)

Parameter	Value
Product Name	ADGRB3 Rabbit pAb
Remarks/Alias	BAI3; ADGRB3
Species	Human
GeneID (Human)	577
GeneID	577
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-280 of human ADGRB3 (NP_001695.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60242
Protein Weight	172kDa
Shipping	Ice bag