

Anti-Adh1B Rabbit Polyclonal Antibody For Wb, Elisa - P00325

Item Number: CM0018185



Introduction

Rabbit polyclonal antibody targeting human ADH1B (all-trans-retinol dehydrogenase), validated for WB and ELISA with rat cross-reactivity. Suitable for retinoid metabolism studies, alcohol metabolism research, and related applications.

[Learn More](#)

Parameter	Value
Product Name	ADH1B Rabbit pAb
Remarks/Alias	ADH2; HEL-S-117; ADH1B
Species	Human
GeneID (Human)	125
GeneID	125
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 206-375 of human ADH1B (NP_000659.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P00325
Protein Weight	40kDa
Shipping	Ice bag