

Anti-Adh5/Gsnor Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P11766

Item Number: CM0023703



Introduction

High-specificity rabbit polyclonal antibody against ADH5/GSNOR (Alcohol dehydrogenase class-3), validated for WB, IF/ICC, and ELISA in human, mouse, and rat. Ideal for formaldehyde detoxification and S-nitrosylation studies.

[Learn More](#)

Parameter	Value
Product Name	ADH5/GSNOR Rabbit pAb
Remarks/Alias	FDH; ADHX; ADH-3; AMEDS; BMFS7; FALDH; GSNOR; GSH-FDH; HEL-S-60p; ADH5/GSNOR
Species	Human
Gene ID (Human)	128
Gene ID	128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-374 of human ADH5/GSNOR (NP_000662.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11766
Protein Weight	40 kDa
Shipping	Ice bag