

Anti-Adh7 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P40394

Item Number: CM0027800



Introduction

Rabbit polyclonal antibody against human ADH7 (retinol dehydrogenase). Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. Suitable for studying retinoid metabolism and aldehyde detoxification.

[Learn More](#)

Parameter	Value
Product Name	ADH7 Rabbit pAb
Remarks/Alias	ADH4; ADH7
Species	Human
GeneID (Human)	131
GeneID	131
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-386 of human ADH7 (NP_000664.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40394
Protein Weight	41kDa
Shipping	Ice bag