

Anti-Aes Rabbit Polyclonal Antibody For Wb, Elisa - Q08117

Item Number: CM0026768

Introduction

Rabbit polyclonal antibody against AES/TLE5, a transcriptional corepressor in the Wnt pathway. Suitable for Western blot and ELISA, reactive in human and mouse. UniProt Q08117.

[Learn More](#)

Parameter	Value
Product Name	AES Rabbit pAb
Remarks/Alias	AES; GRG; ESP1; GRG5; AES-1; AES-2; Grg-5
Species	Human
Gene ID	166 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human AES (NP_945320.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q08117
Protein Weight	22kDa
Shipping	Ice bag