


Anti-Agrin Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O00468

Item Number: CM0023168



Introduction

Rabbit polyclonal antibody against human Agrin validated for Western blot, IF/ICC, and ELISA. Detects human, mouse, and rat Agrin. Ideal for NMJ, cartilage, and synapse research.

[Learn More](#)

Parameter	Value
Product Name	AGRN Rabbit pAb
Remarks/Alias	CMS8; AGRIN; CMSPPD; AGRN
Species	Human
GeneID (human)	375790
GeneID	375790
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-350 of human AGRN (NP_940978.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00468
Protein Weight	217kDa
Shipping	Ice bag