

Anti-Agxt Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P21549

Item Number: CM0025359



Introduction

Rabbit polyclonal antibody against human AGXT (alanine--glyoxylate aminotransferase), validated for WB, IHC-P, ELISA. Cross-reacts with human and mouse. Targets a 43kDa peroxisomal protein involved in glyoxylate detoxification.

[Learn More](#)

Parameter	Value
Product Name	AGXT Rabbit pAb
Remarks/Alias	AGT; PH1; SPT; AGT1; SPAT; TLH6; AGXT1; Ser-PyrAT; AGXT
Species	Human
GeneID (Human)	189
GeneID	189
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P21549
Protein Weight	43kDa
Shipping	Ice bag