

Anti-Ahr Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P35869

Item Number: CM0014386



Introduction

Rabbit polyclonal antibody against human Aryl hydrocarbon receptor (AHR). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Useful for xenobiotic metabolism, immunity, and circadian rhythm studies.

[Learn More](#)

Parameter	Value
Product Name	AHR Rabbit pAb
Remarks/Alias	RP85; bHLHe76; AHR
Species	Human
GeneID (Human)	196
GeneID	196
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35869
Protein Weight	96kDa
Shipping	Ice bag