


Anti-Aire Polyclonal Antibody For Wb, If/Icc, Elisa - O43918

Item Number: CM0020375



Introduction

AIRE rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets autoimmune regulator (APECED protein) involved in thymic self-tolerance. SWISS: O43918.

[Learn More](#)

Parameter	Specification
Product Name	AIRE Rabbit pAb
Remarks/Alias	APS1; APSI; PGA1; AIRE1; APECED; AIRE
Species	Human
Gene ID (Human)	326
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 30-130 of human AIRE (NP_000374.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43918
Protein Weight	58kDa
Shipping	Ice bag