

Anti-Ak7 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96M32

Item Number: CM0015163

Introduction

Rabbit polyclonal antibody against human AK7 (Adenylate kinase 7) for WB, IF/ICC, ELISA. Validated in human, mouse, and rat. Useful for studying ciliary function and nucleotide metabolism.

[Learn More](#)

Parameter	Value
Product Name	AK7 Rabbit pAb
Remarks/Alias	AK 7; FAP75; CFAP75; SPGF27; AK7
Species	Human
Gene ID (Human)	122481
Gene ID	122481
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96M32
Protein Weight	83kDa
Shipping	Ice bag