

Anti-Akap8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43823

Item Number: CM0031066

Introduction

Rabbit polyclonal antibody targeting human AKAP8 (AKAP95) validated for WB, IF/ICC, ELISA; cross-reacts with mouse. Ideal for chromatin condensation, cell cycle regulation, and PKA signaling research.

[Learn More](#)

Parameter	Value
Product Name	AKAP8 Rabbit pAb
Remarks/Alias	AKAP-8; AKAP95; AKAP 95; AKAP-95; AKAP8
Species	Human
GeneID (Human)	10270
GeneID	10270
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-692 of human AKAP8 (NP_005849.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43823
Protein Weight	76kDa
Shipping	Ice bag