

Anti-Akr1A1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P14550

Item Number: CM0017040



Introduction

Rabbit polyclonal antibody against human AKR1A1, validated for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for aldehyde detoxification and drug metabolism research.

[Learn More](#)

Parameter	Value
Product Name	AKR1A1 Rabbit pAb
Remarks/Alias	ALR, ARM, DD3, ALDR1, HEL-S-6, AKR1A1
Species	Human
Gene ID (Human)	10327
Gene ID	10327
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human AKR1A1 (NP_006057.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14550
Protein Weight	37kDa
Shipping	Ice bag