

Anti-Akt1 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P31749

Item Number: CM0016417



Introduction

High-quality AKT1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, IHC-P, and ELISA applications. Detects human, mouse, and rat AKT1 (~56 kDa). Ideal for studies of PI3K/AKT signaling, cell survival, and cancer biology.

[Learn More](#)

Parameter	Description
Product Name	AKT1 Rabbit pAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; AKT1
Species	Human
GeneID (Human)	207
GeneID	207
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 350-480 of human AKT1 (NP_005154.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749
Protein Weight	56 kDa
Shipping	Ice bag