


Anti-Akt2 Polyclonal Antibody For Wb, Elisa - P31751

Item Number: CM0021831



Introduction

Rabbit polyclonal antibody targeting human AKT2 (PKB beta). Validated for WB and ELISA, with cross-reactivity to human and mouse. Suitable for studying insulin signaling, glucose transport, apoptosis, and cancer pathways.

[Learn More](#)

Parameter	Value
Product Name	AKT2 Rabbit pAb
Remarks/Alias	PKBB; PRKBB; HIHGHH; PKBBETA; RAC-BETA; AKT2
Species	Human
GeneID	208
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 381-481 of human AKT2 (NP_001617.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P31751
Protein Weight	56kDa
Shipping	Ice bag