

Anti-Aldh18A1 Rabbit Polyclonal Antibody For Wb, Elisa - P54886

Item Number: CM0016225



Introduction

High-quality ALDH18A1 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, rat. Ideal for studying proline biosynthesis and mitochondrial metabolism.

[Learn More](#)

Parameter	Value
Product Name	ALDH18A1 Rabbit pAb
Remarks/Alias	GSAS; P5CS; PYCS; SPG9; ADCL3; SPG9A; SPG9B; ARCL3A; ALDH18A1
Species	Human
GeneID (Human)	5832
GeneID	5832
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human ALDH18A1 (NP_002851.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54886
Protein Weight	87kDa
Shipping	Ice bag