


## Anti-Aldh2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P05091

Item Number: CM0015125



### Introduction

Rabbit polyclonal antibody against human ALDH2, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying mitochondrial aldehyde dehydrogenase in formaldehyde clearance and related diseases.

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ALDH2 Rabbit pAb                                                                                                                       |
| Remarks/Alias    | ALDM; ALDH1; ALDH-E2; ALDH2                                                                                                            |
| Species          | Human                                                                                                                                  |
| GeneID (human)   | 217                                                                                                                                    |
| GeneID           | 217                                                                                                                                    |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 321-480 of human ALDH2 (NP_000681.2) |
| Source           | Rabbit                                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                                  |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                               |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                      |
| SWISS            | P05091                                                                                                                                 |
| Protein Weight   | 56kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                |