


Anti-Aldh7A1 Monoclonal Antibody For Wb, Elisa - P49419

Item Number: CM0002033



Introduction

Anti-ALDH7A1 Rabbit Monoclonal Antibody validated for WB and ELISA. Detects human and mouse ALDH7A1, a multifunctional aldehyde dehydrogenase involved in oxidative stress defense, lysine catabolism, and ferroptosis inhibition. Ideal for cancer and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ALDH7A1 Rabbit mAb
Remarks/Alias	EPD; PDE; ATQ1; ALDH7A1
Species	Human
GeneID (human)	501
GeneID	501
Immunogen	Synthetic peptide. A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ALDH7A1 (P49419).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49419
Protein Weight	58kDa
Shipping	Ice bag