


Anti-Aldh7A1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49419

Item Number: CM0028906



Introduction

High-quality rabbit polyclonal antibody targeting human ALDH7A1, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying aldehyde metabolism, ferroptosis, and epilepsy.

[Learn More](#)

Parameter	Specification
Product Name	ALDH7A1 Rabbit pAb
Remarks/Alias	EPD; PDE; ATQ1; ALDH7A1
Species	Human
Gene ID	501
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 430-539 of human ALDH7A1 (NP_001173.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49419
Protein Weight	58kDa
Shipping	Ice bag