


Anti-Alkaline Phosphatase (Alpl) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P05186

Item Number: CM0017511



Introduction

Rabbit polyclonal antibody targeting human Alkaline Phosphatase (ALPL), validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for skeletal mineralization and thermogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Alkaline Phosphatase (ALPL) Rabbit pAb
Remarks/Alias	HOPS; HPPA; HPPC; HPPI; HPPO; TNAP; TNALP; APTNAP; TNSALP; AP-TNAP; TNS-ALP; Alkaline Phosphatase (ALPL)
Species	Human
Gene ID (Human)	249
Gene ID	249
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 120-260 of human Alkaline Phosphatase (ALPL) (NP_000469.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05186
Protein Weight	57kDa
Shipping	Ice bag