

## Anti-Alkbh1 Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - Q13686

Item Number: CM0008567



### Introduction

ALKBH1 Rabbit mAb (Q13686) validated for WB, IP, IHC-P, ELISA. Recognizes human, mouse, rat ALKBH1. Suitable for tRNA/DNA demethylation, DNA repair, and mitochondrial translation studies.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | ALKBH1 Rabbit mAb                                                                                                     |
| Remarks/Alias    | ABH; ABH1; alkB; hABH; ALKBH; ALKBH1                                                                                  |
| Species          | Human                                                                                                                 |
| Gene ID (Human)  | 8846                                                                                                                  |
| Gene ID          | 8846                                                                                                                  |
| Immunogen        | Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ALKBH1 (Q13686). |
| Source           | Rabbit                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                 |
| Application      | WB, IP, IHC-P, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                     |
| SWISS            | Q13686                                                                                                                |
| Protein Weight   | 44 kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                               |