

Anti-Alkbh1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q13686

Item Number: CM0013869



Introduction

ALKBH1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, IP, ELISA on human and mouse samples. Targets the nucleic acid dioxygenase ALKBH1, a key enzyme in tRNA demethylation and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	ALKBH1 Rabbit pAb
Remarks/Alias	ABH; ABH1; alkB; hABH; ALKBH; ALKBH1
Species	Human
Gene ID (Human)	8846
Gene ID	8846
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 110-389 of human ALKBH1 (NP_006011.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13686
Protein Weight	44kDa
Shipping	Ice bag