


Anti-Alox12 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P18054

Item Number: CM0015018



Introduction

Rabbit polyclonal antibody against human ALOX12 (12-lipoxygenase). Validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Detects a ~76 kDa protein. Suitable for platelet activation, inflammation, and tumor biology research.

[Learn More](#)

Parameter	Value
Product Name	ALOX12 Rabbit pAb
Remarks/Alias	LOG12; 12-LOX; 12S-LOX; ALOX12
Species	Human
GeneID (Human)	239
GeneID	239
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 483-663 of human ALOX12 (NP_000688.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18054
Protein Weight	76 kDa
Shipping	Ice bag