


## Anti-Alpha-1 Antichymotrypsin (Serpina3) Rabbit Monoclonal Antibody For Wb, Elisa - P01011

Item Number: CM0002051



### Introduction

Rabbit monoclonal antibody against human Alpha-1 Antichymotrypsin (SERPINA3, ACT). Validated for WB and ELISA; targets residues 1-423. Ideal for studies of acute-phase response and amyloid plaque pathology. Recombinant immunogen.

[Learn More](#)

| Parameter        | Value                                                        |
|------------------|--------------------------------------------------------------|
| Product Name     | Alpha-1 Antichymotrypsin (SERPINA3) Rabbit mAb               |
| Remarks/Alias    | ACT; AACT; GIG24; GIG25; Alpha-1 Antichymotrypsin (SERPINA3) |
| Species          | Human                                                        |
| Gene ID (Human)  | 12                                                           |
| Gene ID          | 12                                                           |
| Immunogen        | Recombinant protein                                          |
| Source           | Rabbit                                                       |
| Category         | Monoclonal Antibodies                                        |
| Application      | WB, ELISA                                                    |
| Cross Reactivity | Human                                                        |
| SWISS            | P01011                                                       |
| Protein Weight   | 48kDa                                                        |
| Shipping         | Ice bag                                                      |