


Anti-Alpha-1 Antichymotrypsin (Serpina3) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P01011

Item Number: CM0017500



Introduction

Rabbit polyclonal antibody against human Alpha-1-antichymotrypsin (SERPINA3), validated for WB, IHC-P, ELISA. Inhibits neutrophil cathepsin G and mast cell chymase, implicated in acute phase response and Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	Alpha-1 Antichymotrypsin (SERPINA3) Rabbit pAb
Remarks/Alias	ACT; AACT; GIG24; GIG25; Alpha-1 Antichymotrypsin (SERPINA3)
Species	Human
GeneID (Human)	12
GeneID	12
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-423 of human Alpha-1 Antichymotrypsin (Alpha-1 Antichymotrypsin (SERPINA3)) (NP_001076.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01011
Protein Weight	48kDa
Shipping	Ice bag