

Anti-Alpha-2-Macroglobulin (A2M) Rabbit Polyclonal Antibody For Wb, Elisa - P01023

Item Number: CM0014902



Introduction

Alpha-2-Macroglobulin (A2M) Rabbit Polyclonal Antibody validated for WB and ELISA; detects human, mouse, rat A2M. Used in protease inhibition research. UniProt P01023.

[Learn More](#)

Parameter	Value
Product Name	Alpha-2-Macroglobulin (A2M) Rabbit pAb
Remarks/Alias	A2MD; CPAMD5; FWP007; S863-7; Alpha-2-Macroglobulin (A2M)
Species	Human
GeneID (Human)	2
GeneID	2
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1155-1474 of human Alpha-2-Macroglobulin (A2M) (NP_000005.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01023
Protein Weight	163kDa
Shipping	Ice bag