

Anti-Alpha-Fetoprotein (Afp) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P02771

Item Number: CM0005840



Introduction

Alpha-Fetoprotein (AFP) Rabbit Monoclonal Antibody, produced in rabbit, recognizes human AFP. Suitable for WB, IF/ICC, ELISA applications. Targets a 69kDa secreted biomarker synthesized by fetal liver.

[Learn More](#)

Parameter	Value
Product Name	Alpha-Fetoprotein (AFP) Rabbit mAb
Remarks/Alias	AFPD; FETA; HPAFP; Alpha-Fetoprotein (AFP)
Species	Human
GeneID	174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 395-609 of human Alpha-Fetoprotein (AFP) (NP_001125.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P02771
Protein Weight	69kDa
Shipping	Ice bag