


## Anti-Ambra1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9C0C7

Item Number: CM0018712



### Introduction

Rabbit polyclonal antibody targeting human AMBRA1 (DCAF3/WDR94), validated for WB and ELISA, with mouse cross-reactivity. Detects the 143 kDa autophagy and cell cycle regulator. Ideal for IVD research and assay development.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Ambra1 Rabbit pAb                                                                                                    |
| Remarks/Alias    | DCAF3; WDR94; Ambra1                                                                                                 |
| Species          | Human                                                                                                                |
| GeneID (Human)   | 55626                                                                                                                |
| GeneID           | 55626                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 961-1208 of human Ambra1 (NP_060219.2) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Mouse                                                                                                                |
| SWISS            | Q9C0C7                                                                                                               |
| Protein Weight   | 143kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                              |