


Anti-Angptl3 Polyclonal Antibody For Wb, Elisa - Q9Y5C1

Item Number: CM0014448



Introduction

Rabbit polyclonal antibody targeting human ANGPTL3, recommended for Western blotting and ELISA, cross-reactive with mouse. Suitable for lipid metabolism and angiogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	ANGPTL3 Rabbit pAb
Remarks/Alias	ANL3, ANG-5, FHBL2, ANGPT5, ANGPTL3
Species	Human
GeneID	27329 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 231-460 of human ANGPTL3 (NP_055310.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9Y5C1
Protein Weight	54kDa
Shipping	Ice bag