

Anti-Ap2A1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O95782

Item Number: CM0011283

Introduction

Recombinant rabbit monoclonal antibody targeting human AP2A1 (SWISS O95782) for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for clathrin-mediated endocytosis and membrane trafficking research.

[Learn More](#)

Parameter	Value
Product Name	AP2A1 Rabbit mAb
Remarks/Alias	ADTAA; CLAPA1; AP2-ALPHA; AP2A1
Species	Human
GeneID (Human)	160
GeneID	160
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human AP2A1 (O95782).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95782
Protein Weight	108kDa
Shipping	Ice bag