


Anti-Ap2S1 Rabbit Polyclonal Antibody For Wb And Elisa - P53680

Item Number: CM0015431



Introduction

Rabbit polyclonal antibody against AP2S1 (AP-2 complex subunit sigma) validated for WB and ELISA. Recognizes human, mouse, and rat AP2S1. Immunogen is a recombinant protein covering the N-terminal region. Essential for studying clathrin-mediated endocytosis and protein trafficking.

[Learn More](#)

Parameter	Specification
Product Name	AP2S1 Rabbit pAb
Remarks/Alias	AP17; FBH3; HHC3; FBHOK; CLAPS2; AP2S1
Species	Human
Gene ID (Human)	1175
Gene ID	1175
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human AP2S1 (NP_004060.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53680
Protein Weight	17kDa
Shipping	Ice bag