

Anti-AP3S2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P59780

Item Number: CM0021560

Introduction

Anti-AP3S2 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat AP3S2, a subunit of the AP-3 complex involved in Golgi-to-lysosome vesicle trafficking. For cellular transport and lysosomal studies.

[Learn More](#)

Parameter	Value
Product Name	AP3S2 Rabbit pAb
Remarks/Alias	AP3S3; sigma3b; AP3S2
Species	Human
Gene ID (Human)	10239
Gene ID	10239
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human AP3S2 (NP_005820.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P59780
Protein Weight	22kDa
Shipping	Ice bag