


Anti-Apoc2 Rabbit Polyclonal Antibody For Wb, Elisa - P02655

Item Number: CM0020026



Introduction

High-quality rabbit polyclonal antibody targeting human Apolipoprotein C-II (APOC2). Validated for WB and ELISA applications with cross-reactivity to human and rat. Ideal for lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	APOC2 Rabbit pAb
Remarks/Alias	APO-CII; APOC-II; APOC2
Species	Human
Gene ID (Human)	344
Gene ID	344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-101 of human APOC2 (NP_000474.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P02655
Protein Weight	11kDa
Shipping	Ice bag