

Anti-Apoe Monoclonal Antibody For Wb, If/Icc, Elisa - P02649

Item Number: CM0009084



Introduction

High-quality rabbit monoclonal anti-APOE antibody validated for WB, IF/ICC, ELISA. Recognizes human APOE protein (36 kDa) involved in lipid transport and Alzheimer's disease research. Ideal for IVD and research.

[Learn More](#)

Parameter	Value
Product Name	APOE Rabbit mAb
Remarks/Alias	AD2, LPG, APO-E, ApoE4, LDLCQ5
Species	Human
GeneID (Human)	348
GeneID	348
Immunogen	Recombinant protein corresponding to amino acids 19-317 of human APOE (NP_000032.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P02649
Protein Weight	36kDa
Shipping	Ice bag