

Anti-Apoer2/Lrp8 Rabbit Monoclonal Antibody For Wb, Elisa - Q14114

Item Number: CM0004609



Introduction

High-affinity rabbit monoclonal antibody targeting human ApoER2/LRP8, suitable for Western blotting and ELISA. Detects human, mouse, and rat LRP8. Ideal for Reelin signaling, neuronal development, and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	ApoER2/LRP8 Rabbit mAb
Remarks/Alias	MC11; LRP-8; APOER2; HSZ75190; ApoER2/LRP8
Species	Human
Gene ID (Human)	7804
Gene ID	7804
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 864-963 of human ApoER2/LRP8 (Q14114).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14114
Protein Weight	106kDa
Shipping	Ice bag