

Anti-APOL1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O14791

Item Number: CM0012273



Introduction

Rabbit monoclonal antibody against human APOL1, suitable for WB, IF/ICC, ELISA. APOL1 is involved in lipid transport and trypanosome lysis, ideal for metabolism and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	APOL1 Rabbit mAb
Remarks/Alias	APOL; APO-L; FSGS4; APOL-I; APOL1
Species	Human
GeneID (Human)	8542
GeneID	8542
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 303-398 of human APOL1 (O14791)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O14791
Protein Weight	44kDa
Shipping	Ice bag