

Anti-Aqp2 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P41181

Item Number: CM0021623



Introduction

AQP2 Rabbit pAb is a polyclonal antibody for WB, IF-P, IHC-P, ELISA; detects human AQP2 with cross-reactivity to mouse, rat. Suitable for renal water transport and diabetes insipidus studies.

[Learn More](#)

Parameter	Value
Product Name	AQP2 Rabbit pAb
Remarks/Alias	NDI2; AQP-CD; WCH-CD; AQP2
Species	Human
GeneID (Human)	359
GeneID	359
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 147-271 of human AQP2 (NP_000477.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41181
Protein Weight	29kDa
Shipping	Ice bag