

Anti-Aquaporin-8 (Aqp8) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O94778

Item Number: CM0028816

Introduction

Rabbit polyclonal anti-AQP8 antibody validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for studying aquaporin-8 water and hydrogen peroxide transport, bile secretion, and mitochondrial function. Predicted molecular weight ~27 kDa.

[Learn More](#)

Parameter	Value
Product Name	AQP8 Rabbit pAb
Remarks/Alias	AQP-8; AQP8
Species	Human
Gene ID (Human)	343
Gene ID	343
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human AQP8 (NP_001160.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94778
Protein Weight	27kDa
Shipping	Ice bag